



Email : sales@laboratoryinstrumentindia.com

Product Name :

Proshield 10 Basic Disposable Lab Coats

Product Code :

LBNY-00017-11000357



Description :

Disposable lab coats that are optimized for comfort, softness and breathability
Spun bond melt blown spun bond (SMS) material
Designed for workers in a range of industries for nonhazardous dry particle applications.
Snap closure (six plus one adjustable)
Serged seams
Three pockets (one left chest pencil, two lower front)
Set sleeves
Knit cuffs.
Knit collar.

Technical Specification :

Proshield 10 Basic Disposable Lab Coats



Laboratory Instrument India