



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hycar Aprons

Product Code :
LBNY-00017-11000359



Description :

Hycar Aprons keep workers clean and protected from harmful contact wherever oil, grease and chemicals are present

Construction is traditionally tougher and more cut, abrasion and puncture resistant than neoprene but is still flexible for the toughest jobs

Coated poly/cotton reinforcement is easy to clean

Sealed grommets allow for size adjustment without tear-out

Tough nitrile blend

Ideal for grease and animal fats

Excellent in cold conditions

Bib styles with adjustable ties ensure optimal fit for a wide range of users

Applications:

Storing and transportation

Equipment maintenance and repair

Plant disinfection and sanitation

Packing of fillets and ground meat

Canning; Bottling

Transferring liquids.

Technical Specification :

Hycar Aprons



Laboratory Instrument India