



Email : sales@laboratoryinstrumentindia.com

Product Name :
Endurosaf Aprons

Product Code :
LBNY-00017-11000360



Description :

High performance, adaptable protection for a wide variety of jobs and are resistant to abrasion, chemicals, oils and fats

Soft, comfortable, easy to clean and remain flexible in cold temperatures.

Strong, tough, 8 to 9 mil die cut polyurethane material offers a unique combination of characteristics to outperform neoprene, nitrile, vinyl and other films

Sealed ties and adjustable neck straps provide a durable, consistent fit for a wide variety of users.

Technical Specification :

Endurosaf Aprons



Laboratory Instrument India