



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Endurosaf Aprons

**Product Code :**  
LBNY-00017-11000360



**Description :**

High performance, adaptable protection for a wide variety of jobs and are resistant to abrasion, chemicals, oils and fats

Soft, comfortable, easy to clean and remain flexible in cold temperatures.

Strong, tough, 8 to 9 mil die cut polyurethane material offers a unique combination of characteristics to outperform neoprene, nitrile, vinyl and other films

Sealed ties and adjustable neck straps provide a durable, consistent fit for a wide variety of users.

**Technical Specification :**

Endurosaf Aprons



**Laboratory Instrument India**