



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cut Resistant Sleeve

**Product Code :**  
LBNY-00017-11000361



**Description :**

Made from 7 gauge knitted heavyweight Spectra it offers high level cut protection for the shoulder  
Made from 7 gauge knitted heavyweight Spectra™ this sleeve offers high level cut protection up to the shoulder  
Spectra yarns are wrapped in polyester which offers better comfort and flexibility  
Extremely cut resistant, flexible and comfortable  
An attached adjustable strap is included to hold the sleeve in place  
Can be bleached while laundering  
Uses:  
Automotive, Food processing, Metal fabrication.

**Technical Specification :**

Cut Resistant Sleeve



**Laboratory Instrument India**