



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

10000FR Encapsulated Level A Suit with MSA Quick Fill

**Product Code :**

LBNY-00017-11000366



**Description :**

10000 FR fabrics are tested by a third-party laboratory to assess its chemical permeation resistance. Provides at least 30 minutes of protection against over 285 chemical challenges  
10000 FR is a top choice for chemical protection in emergency situations.  
Certified to NFPA 1991 in a single-skin construction, 10000 FR garments are composed of a patented, heavy-duty aluminized fabric made of a blend, laminated to a multi-layer barrier film.  
10000 FR is high-visibility silver  
Ideally suited for industrial and HAZMAT situations  
A limited-use "single skin" garment offers broad chemical holdout.

**Technical Specification :**

10000FR Encapsulated Level A Suit with MSA Quick Fill



**Laboratory Instrument India**