



Email : sales@laboratoryinstrumentindia.com

Product Name :

10000FR Encapsulated Level A Suit with MSA Quick Fill

Product Code :

LBNY-00017-11000366



Description :

10000 FR fabrics are tested by a third-party laboratory to assess its chemical permeation resistance. Provides at least 30 minutes of protection against over 285 chemical challenges
10000 FR is a top choice for chemical protection in emergency situations.
Certified to NFPA 1991 in a single-skin construction, 10000 FR garments are composed of a patented, heavy-duty aluminized fabric made of a blend, laminated to a multi-layer barrier film.
10000 FR is high-visibility silver
Ideally suited for industrial and HAZMAT situations
A limited-use "single skin" garment offers broad chemical holdout.

Technical Specification :

10000FR Encapsulated Level A Suit with MSA Quick Fill



Laboratory Instrument India