



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
10000FR Encapsulated Level A Suit with Hansen Fitting

**Product Code :**  
LBNY-00017-11000370



**Description :**

Provides at least 30 minutes of protection against over 285 chemical challenges  
Ideally suited for industrial  
10000 FR is high-visibility silver  
10000 FR is a top choice for chemical protection in emergency situations.  
10000 FR fabrics are tested by a third-party laboratory to assess its chemical permeation resistance.  
A limited-use “single skin” garment offers broad chemical holdout.

**Technical Specification :**

10000FR Encapsulated Level A Suit with Hansen Fitting



**Laboratory Instrument India**