



Email : sales@laboratoryinstrumentindia.com

Product Name :
10000FR Encapsulated Level A Suit with Hansen Fitting

Product Code :
LBNY-00017-11000370



Description :

Provides at least 30 minutes of protection against over 285 chemical challenges
Ideally suited for industrial
10000 FR is high-visibility silver
10000 FR is a top choice for chemical protection in emergency situations.
10000 FR fabrics are tested by a third-party laboratory to assess its chemical permeation resistance.
A limited-use "single skin" garment offers broad chemical holdout.

Technical Specification :

10000FR Encapsulated Level A Suit with Hansen Fitting



Laboratory Instrument India