



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Denim Shop Aprons, Bib Pocket and One Side Pocket

**Product Code :**  
LBNY-00017-11000375



**Description :**

Preshrunk 100% cotton, 14 oz. denim minimizes shrinkage  
Non pre-laced polyester ties for neck and waist  
Available in choice of several pocket styles  
With metal grommets  
Measure 36L x 28 in. W (91 x 71cm).  
Color: Blue.

**Technical Specification :**

Denim Shop Aprons, Bib Pocket and One Side Pocket



**Laboratory Instrument India**