



Email : sales@laboratoryinstrumentindia.com

Product Name :
Denim Shop Aprons, Bib Pocket and One Side Pocket

Product Code :
LBNY-00017-11000375



Description :

Preshrunk 100% cotton, 14 oz. denim minimizes shrinkage
Non pre-laced polyester ties for neck and waist
Available in choice of several pocket styles
With metal grommets
Measure 36L x 28 in. W (91 x 71cm).
Color: Blue.

Technical Specification :

Denim Shop Aprons, Bib Pocket and One Side Pocket



Laboratory Instrument India