



Email : sales@laboratoryinstrumentindia.com

Product Name :
2000 Series 125 Coveralls

Product Code :
LBNY-00017-11000377



Description :

Flexible and durable light liquid splash protection
Elastic opening for tighter fit at wrist
Elastic opening for tighter fit at ankle
QC is polyethylene-coated fabric provides at least 30 minutes of protection against 42 chemical challenges
Laydown collar
10 mil thickness
Yellow for high visibility.

Technical Specification :

2000 Series 125 Coveralls



Laboratory Instrument India