



Email : sales@laboratoryinstrumentindia.com

Product Name :
CPF 3 Aprons

Product Code :
LBNY-00017-11000378



Description :

For rigorous activities and rugged situations involving liquid splash
CPF 3 provides at least 30 minutes of protection against 149 chemical challenges
Taped seams and bound snap neck
CPF 3 is multi-layer barrier film laminated to durable 2.0 oz./yd² polypropylene fabric
Strong and durable 44 in. long apron with broad chemical barrier
Apron with attached long sleeves, waist ties and elastic wrists
Tan color for discretionary purposes.

Technical Specification :

CPF 3 Aprons



Laboratory Instrument India