



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
CPF 3 Aprons

**Product Code :**  
LBNY-00017-11000378



**Description :**

For rigorous activities and rugged situations involving liquid splash  
CPF 3 provides at least 30 minutes of protection against 149 chemical challenges  
Taped seams and bound snap neck  
CPF 3 is multi-layer barrier film laminated to durable 2.0 oz./yd<sup>2</sup> polypropylene fabric  
Strong and durable 44 in. long apron with broad chemical barrier  
Apron with attached long sleeves, waist ties and elastic wrists  
Tan color for discretionary purposes.

**Technical Specification :**

CPF 3 Aprons



**Laboratory Instrument India**