



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
50 Coveralls with Elastic Wrist and Ankles

**Product Code :**  
LBNY-00017-11000383



**Description :**

Made with a microporous film that is laminated to a spund bonded polypropylene substrate  
Front zipper closure  
Made from microporous film, Pro Shield 50 garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.  
Elastic ankles  
Serged seams  
Elastic wrist  
Collar  
White  
Thickness: 10 mils.

**Technical Specification :**

50 Coveralls with Elastic Wrist and Ankles



**Laboratory Instrument India**

