



Email : sales@laboratoryinstrumentindia.com

Product Name :
50 Coveralls with Elastic Wrist and Ankles

Product Code :
LBNY-00017-11000383



Description :

Made with a microporous film that is laminated to a spund bonded polypropylene substrate
Front zipper closure
Made from microporous film, Pro Shield 50 garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.
Elastic ankles
Serged seams
Elastic wrist
Collar
White
Thickness: 10 mils.

Technical Specification :

50 Coveralls with Elastic Wrist and Ankles



Laboratory Instrument India

