



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
White Shirt

**Product Code :**  
LBNY-00017-11000388



**Description :**

Shirt design with front snap closures for easy donning and doffing

Shirt with laydown collar and open wrists

Front snap closure

Serged seam and shirt extends to hip

Garments are composed of flash spun high density polyethylene which creates a unique, nonwoven material. It provides an ideal balance of protection, durability and comfort of any limited use fabric technology.

The fabric offers an inherent barrier against particles (down to 1.0 micron in size).

Applications

Lead and asbestos abatement/remediation, general maintenance/operations, spray painting, general clean-up.

**Technical Specification :**

White Shirt



**Laboratory Instrument India**