



Email : sales@laboratoryinstrumentindia.com

Product Name :
NexGen Apron

Product Code :
LBNY-00017-11000389



Description :

Maximum protection at affordable cost
Bib-style apron provides coverage of front and sides
With bound neck and ties
NexGen is liquid resistant and provides a true aerosol barrier, with >99% holdout of particles down to 0.35µm.
Proprietary microporous technology gives aprons effective protective properties at an economical price. Perfect for applications involving nontoxic liquid or spray hazards.
Apron
Length: 36 in. (91cm)
White.

Technical Specification :

NexGen Apron



Laboratory Instrument India