



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
NexGen Apron

**Product Code :**  
LBNY-00017-11000389



**Description :**

Maximum protection at affordable cost  
Bib-style apron provides coverage of front and sides  
With bound neck and ties  
NexGen is liquid resistant and provides a true aerosol barrier, with >99% holdout of particles down to 0.35µm.  
Proprietary microporous technology gives aprons effective protective properties at an economical price. Perfect for applications involving nontoxic liquid or spray hazards.  
Apron  
Length: 36 in. (91cm)  
White.

**Technical Specification :**

NexGen Apron



**Laboratory Instrument India**