



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coveralls, Clean-Processed and Sterile

Product Code :
LBNY-00017-11000393



Description :

Lightweight, breathable, and clean processed coveralls for use in ISO 4-5 Class 10-100 critical environments. Material and bound seams offer an excellent barrier while maintaining comfort and durability
Gamma sterilized and clean processed
Spun polyethylene creates barrier to particles, micro-organisms and non-hazardous liquids
For use in controlled environments with high standards for particle and microbiological contamination control (pharmaceutical, biotech, medical device and electronic settings).
Lowest level particle shedding
Bound neck and seams
Elastic wrists and ankles
Should not be worn around open flame, heat source, sparks or ignition sources.
Dolman sleeves
Zipper closure

Technical Specification :

Coveralls, Clean-Processed and Sterile



Laboratory Instrument India

