



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Coveralls, Clean-Processed and Sterile

**Product Code :**  
LBNY-00017-11000393



**Description :**

Lightweight, breathable, and clean processed coveralls for use in ISO 4-5 Class 10-100 critical environments. Material and bound seams offer an excellent barrier while maintaining comfort and durability  
Gamma sterilized and clean processed  
Spun polyethylene creates barrier to particles, micro-organisms and non-hazardous liquids  
For use in controlled environments with high standards for particle and microbiological contamination control (pharmaceutical, biotech, medical device and electronic settings).  
Lowest level particle shedding  
Bound neck and seams  
Elastic wrists and ankles  
Should not be worn around open flame, heat source, sparks or ignition sources.  
Dolman sleeves  
Zipper closure

**Technical Specification :**

Coveralls, Clean-Processed and Sterile



**Laboratory Instrument India**

