



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
A8 Zippered Lab Coat

**Product Code :**  
LBNY-00017-11000396



**Description :**

Soft knit cuffs

High mandarin collar

Double bagged.

A8 Zippered Lab Coats are made with a protective 3-layer SMS fabric and are recommended where there is the potential for light fluid contact.

Zippered front closure and no pockets to reduce risk of contamination.

**Technical Specification :**

A8 Zippered Lab Coat



**Laboratory Instrument India**