



Email : sales@laboratoryinstrumentindia.com

Product Name :
10,000 Level A Suits

Product Code :
LBNY-00017-11000400



Description :

Fully Encapsulated Level A Suits are ideal for industrial, HazMat and domestic preparedness situations. High-visibility lime yellow.
Three-layer face shield (PVC 40mil/PTFE 5mil/PVC 20mil)
Gas-tight PVC zipper closure
Double storm flaps with Velcro closure
Double taped seams
Two exhaust valves
Internal adjustment belt
Garments successfully tested against more than 260 challenge chemicals
Attached internal Barrier gloves, outer butyl gloves and sock boots with outer boot flaps
Knee wear pads.

Technical Specification :

10,000 Level A Suits



Laboratory Instrument India

