



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sterile Coveralls with Collar and Thumb Loop

Product Code :
LBNY-00017-11000402



Description :

Disposable sterile coverall with thumb loop for secure hold.
Low linting, anti-static Clean Tough Fabric made
Zip front with flap cover
Sterilized by Gamma irradiation
Anti-Static properties
Low-particulate
Sterility Assurance Level 10⁻⁶.
Elasticated back, cuffs and ankles
Available in sizes ranging from S through 4XL.

Technical Specification :

Sterile Coveralls with Collar and Thumb Loop



Laboratory Instrument India