



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polyethylene Isolation Gown

Product Code :
LBNY-00017-11000409



Description :

Designed to protect healthcare professionals and patients against cross-contamination
Features comfortable fit with thumb loop cuffs
These Polyethylene Isolation Gown are designed to protect health professionals and patients from cross infection. Single-use over-the-head design with thumb loop cuffs.
Made of spunbond polypropylene, latex-free.
Sewn seams provide superior strength
Single-use disposable.

Technical Specification :

Polyethylene Isolation Gown



Laboratory Instrument India