



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Polyethylene Isolation Gown

**Product Code :**  
LBNY-00017-11000409



**Description :**

Designed to protect healthcare professionals and patients against cross-contamination  
Features comfortable fit with thumb loop cuffs  
These Polyethylene Isolation Gown are designed to protect health professionals and patients from cross infection. Single-use over-the-head design with thumb loop cuffs.  
Made of spunbond polypropylene, latex-free.  
Sewn seams provide superior strength  
Single-use disposable.

**Technical Specification :**

Polyethylene Isolation Gown



**Laboratory Instrument India**