



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safe Wear High Performance Lab Coat, Deep Blue

**Product Code :**  
LBNY-00017-11000414



**Description :**

Made from high-quality 45 g/m2 SMS fabric  
Detailed finishing, 3-pocket design and front snap closure  
The Safe Wear High Performance Lab Coat, Deep Blue is premium protective apparel made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.  
Each bag contains 10+2 units, for a total of 12 vs. 10 per package from the competition  
Available in a wide range of styles, colors and sizes  
Comply with recommended OSHA best practices.

**Technical Specification :**

Safe Wear High Performance Lab Coat, Deep Blue



**Laboratory Instrument India**