



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electric Heating Mantle and Electromagnetic Stirrer 1000ml

**Product Code :**  
LBNY-0001-1000013



**Description :**

Electric Heating Mantle and Electromagnetic Stirrer 1000ml

**Technical Specification :**

It is used for liquid stirring and heating in colleges, petrochemical industries, pharmaceutical and environment protection institutions.

Electric Heating Mantle and Electromagnetic Stirrer 1000ml

Digital display, convenient & accurate to control temperature and easy to read out.

With alkali-free glass fiber mantle, corrosion-resistant, heat-resistant and durable.

Fast heating-up without open flame, long lifespan.

Inner hemispherical heater with larger heating area for uniform and efficient heat.

Comes with magnetic stirring function, making the liquid heat evenly under stirring, ideal heating mantle for labs and chemical institutions.



**Laboratory Instrument India**