



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Yellow Coveralls with Hood

**Product Code :**  
LBNY-00017-11000419



**Description :**

These coveralls are designed for applications such as environmental clean-up, sewage purification installations, and composites manufacturing.

Provide an excellent barrier to harmful chemicals and blood borne pathogens.

Yellow Coveralls with Hood are lightweight and relatively strong and durable.

Design: Elastic waist, wrists and ankles, 2-way front zipper with taped storm flap, bound seams, finger loops, 3 piece respirator fit hood

Applications: Environmental cleanup, sewage purification installations, composites manufacturing

Protection: A barrier to numerous inorganic liquid chemicals including acids and bases; achieves the highest classifications for protection from biological for penetration of blood, bodily fluids, and blood-borne pathogens

Comfort: Lightweight yet strong and durable

Sizes: Available in Small through 5X-Large.

**Technical Specification :**

Yellow Coveralls with Hood



**Laboratory Instrument India**

