



Email : sales@laboratoryinstrumentindia.com

Product Name :
Green Coveralls with Hood

Product Code :
LBNY-00017-11000420



Description :

Provide an exceptional barrier against organic and inorganic chemicals as well as biological agents.
Design: Elastic waist, wrists and ankles, double zip closure, double cuffs with knitted inner cuff, ultrasonically welded and taped seams, 2 piece respirator fit hood
Applications: Chemicals handling/transportation, oil-based mud protection, hazardous waste remediation, sewage purification installations, industrial/tank cleaning, HAZMAT emergency response (Level B), pharmaceutical, mining, agriculture
Protection: Permeation tested against over 190 chemicals, including chemical warfare agents; fully encapsulated versions also available
Comfort: Textile like interior improves wearer acceptance
Green Coveralls with Hood incorporate unique design features that enhance wearer comfort and safety.
Sizes: Available in Small through 5X-Large.

Technical Specification :

Green Coveralls with Hood



Laboratory Instrument India

