



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Green Coveralls with Hood

**Product Code :**  
LBNY-00017-11000420



**Description :**

Provide an exceptional barrier against organic and inorganic chemicals as well as biological agents.  
Design: Elastic waist, wrists and ankles, double zip closure, double cuffs with knitted inner cuff, ultrasonically welded and taped seams, 2 piece respirator fit hood  
Applications: Chemicals handling/transportation, oil-based mud protection, hazardous waste remediation, sewage purification installations, industrial/tank cleaning, HAZMAT emergency response (Level B), pharmaceutical, mining, agriculture  
Protection: Permeation tested against over 190 chemicals, including chemical warfare agents; fully encapsulated versions also available  
Comfort: Textile like interior improves wearer acceptance  
Green Coveralls with Hood incorporate unique design features that enhance wearer comfort and safety.  
Sizes: Available in Small through 5X-Large.

**Technical Specification :**

Green Coveralls with Hood



**Laboratory Instrument India**

