



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ammeter

**Product Code :**  
LBNY-0003-1010001



**Description :**

Ammeter

**Technical Specification :**

Square Clear acrylic front cover, fully enclosed economically priced and offer + 2.5% accuracy.  
Moving coil, centre pole movements in clear plastic case.  
Double scaled 0 - 1A, 0 – 5A, with slopping front.



**Laboratory Instrument India**