



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hood with Attached Mask

Product Code :
LBNY-00017-11000425



Description :

Made using a patented flash spinning process, provides an inherent barrier to particles, microorganisms and non-hazardous light liquid splash
Delivers an ideal balance of protection, durability, and comfort.
Hood with attached 7 in. blue mask
Sterile garments are gamma sterilized.
Pleated outer mask is made of polyethylene
Seams and head opening are bound; ties and loops provide a secure fit
Individually packaged in an opaque bag and placed in double transparent poly liners.

Technical Specification :

Hood with Attached Mask



Laboratory Instrument India