



Email : sales@laboratoryinstrumentindia.com

Product Name :
Apron

Product Code :
LBNY-00017-11000429



Description :

Apron is liquid resistant and provides a true aerosol barrier, with >99% holdout of particles down to 0.35µm. Proprietary microporous technology gives aprons effective protective properties at an economical price. Perfect for applications involving nontoxic liquid or spray hazards.

With bound neck and ties

Bib-style apron provides coverage of front and sides

Length: 36 in. (91cm)

White.

Technical Specification :

Apron



Laboratory Instrument India