



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
275 Series Aprons

**Product Code :**  
LBNY-00017-11000436



**Description :**

4000 fabric provides at least 30 minutes of protection against >124 chemical challenges. Typical applications: chemical mixing, hazardous remediation, emergency medical response, paint spraying, and nuclear environments.

Knee-length apron/gown design with attached sleeves Provides enhanced frontal protection

4000 delivers effective protection against a range of chemical environments.

Utilizes film laminated to fabric to provide a rugged and durable fabric that resists abrasion.

Open back for wearer comfort Bound seams around neck opening with snap closure

Ties at waist Elastic opening for tighter fit at wrist 44 in. long for neck-to-knee coverage

Taped seams provide strong chemical resistance against heavy liquid splashes. A sewn seam is covered with a strip of compatible chemical-resistant material through heat sealing.

Available in multiple sizes.

**Technical Specification :**

275 Series Aprons



**Laboratory Instrument India**

