



Email : sales@laboratoryinstrumentindia.com

Product Name :

Safe Wear High Performance Lab Coat, Pretty Pink

Product Code :

LBNY-00017-11000440



Description :

The Safe Wear High Performance Lab Coat, Pretty Pink is premium protective apparel made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style. Each bag contains 10+2 units, for a total of 12 vs. 10 per package from the competition. Fluid-resistant fabric meets AAMI PB70 Level 2 standards. Available in a wide range of styles, colours and sizes. Made from high-quality 45 g/m² SMS fabric. Detailed finishing, 3-pocket design and front snap closure.

Technical Specification :

Safe Wear High Performance Lab Coat, Pretty Pink



Laboratory Instrument India