



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bunsen Burner

**Product Code :**  
LBNY-0003-1010003



**Description :**

Bunsen Burner

**Technical Specification :**

Nickel plated burner tube with rotatable air regulator on enamelled die casted base with a tapered ruffled gas inlet.

Nickel plated burner tube with air regulator, 60 x 10mm (height x diameter) base-height 80mm.



**Laboratory Instrument India**