



Email : sales@laboratoryinstrumentindia.com

Product Name :
1500 PLUS Stitched

Product Code :
LBNY-00017-11000449



Description :

Air and water vapor permeability: Ensures heat stress protection
Particle-filtering, anti-static SMS protective suit, delivering light chemical resistance.
Optimized body fit: Enhanced comfort.
Particle-filtering SMS fabric: Filters 100% of particles over 3 microns(According to particle penetration test).

Technical Specification :

1500 PLUS Stitched



Laboratory Instrument India