



Email : sales@laboratoryinstrumentindia.com

Product Name :
Clean Gown

Product Code :
LBNY-00017-11000454



Description :

Made using a patented flash spinning process
Serged seams have multiple interlocking threads that are sewn around the raw edges of material
Knit cuffs for comfort
Gown provides an inherent barrier to particles, microorganisms and non-hazardous light liquid splash.
Delivers an ideal balance of protection, durability and comfort
One size fits most
Bulk packaged in double transparent poly liners.
Bound neck for lower particle shedding with ties for an adjustable fit
Bound ties at the waist originate at the front and tie in the back

Technical Specification :

Clean Gown



Laboratory Instrument India