



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sleeve with Intercept Technology

Product Code :
LBNY-00017-11000458



Description :

Sleeve with Intercept Technology features a knitted bicep cuff to keep the sleeve in place
Made of black yarn
Silicone-free for friendly paint and finish process
Available in different lengths and widths for optimal fit and comfort to fit.
Color-coding helps for easy sorting after laundry

Applications:

Press shop operations
Glass repair and replacement
Seizing structural and body parts
Body surface checking, marking and touching up
Adjusting systems and screwing
Adjusting panels parts components.
Handling metal sheets and panels
Handling parts with sharp, rough edges
Handling structural and body components
Inspection, selecting, checking parts.

Technical Specification :

Sleeve with Intercept Technology



Laboratory Instrument India