



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sleeve with Intercept Technology

**Product Code :**  
LBNY-00017-11000458



**Description :**

Sleeve with Intercept Technology features a knitted bicep cuff to keep the sleeve in place  
Made of black yarn  
Silicone-free for friendly paint and finish process  
Available in different lengths and widths for optimal fit and comfort to fit.  
Color-coding helps for easy sorting after laundry

**Applications:**

Press shop operations  
Glass repair and replacement  
Seizing structural and body parts  
Body surface checking, marking and touching up  
Adjusting systems and screwing  
Adjusting panels parts components.  
Handling metal sheets and panels  
Handling parts with sharp, rough edges  
Handling structural and body components  
Inspection, selecting, checking parts.

**Technical Specification :**

Sleeve with Intercept Technology



**Laboratory Instrument India**