



Email : sales@laboratoryinstrumentindia.com

Product Name :
Chemical Protective Apron with Sleeves

Product Code :
LBNY-00017-11000461



Description :

Chemical Protective Apron with Sleeves for exceptional chemical protection
Double cuff with knitted inner cuff
Welded and taped seams
Velcro fastening to neck
Tie fastening at the waist

Applications:

Tranferring liquids and solids
Opening and draining pumps valves or lines
Contaminated environments
Chemical handling
Filling, blending and charging of raw materials
Raw material sample collection
Concentration mixing and palletisation
Loading and unloading of process equipment
Assembly and disassembly of large pipes.
Hazmat response.

Technical Specification :

Chemical Protective Apron with Sleeves



Laboratory Instrument India