



Email : sales@laboratoryinstrumentindia.com

Product Name :
50 Coveralls

Product Code :
LBNY-00017-11000465



Description :

Made with a microporous film that is laminated to a spund bonded polypropylene substrate
Made from microporous film, Pro Shield 50 garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.
Front zipper closure
Serged seams
Collar
Open wrist
Open ankles
White
Thickness: 10 mils.

Technical Specification :

50 Coveralls



Laboratory Instrument India

