



Email : sales@laboratoryinstrumentindia.com

Product Name :
182B Coveralls, Clean-Processed

Product Code :
LBNY-00017-11000468



Description :

Made using a patented flash spinning process
Provides an inherent barrier to particles, microorganisms and non-hazardous light liquid splash
Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration
Bound neck for lower particle shedding
Front zipper closure for easy donning and doffing
Clean delivers an ideal balance of protection, durability, and comfort. Clean-processed garments offer lowest level of particle shedding within DuPont product portfolio.
Elastic openings for tighter fit at wrist and ankle
Raglan sleeve design for greater range of motion and comfort.

Technical Specification :

182B Coveralls, Clean-Processed



Laboratory Instrument India

