



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
182B Coveralls, Clean-Processed

**Product Code :**  
LBNY-00017-11000468



**Description :**

Made using a patented flash spinning process  
Provides an inherent barrier to particles, microorganisms and non-hazardous light liquid splash  
Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration  
Bound neck for lower particle shedding  
Front zipper closure for easy donning and doffing  
Clean delivers an ideal balance of protection, durability, and comfort. Clean-processed garments offer lowest level of particle shedding within DuPont product portfolio.  
Elastic openings for tighter fit at wrist and ankle  
Raglan sleeve design for greater range of motion and comfort.

**Technical Specification :**

182B Coveralls, Clean-Processed



**Laboratory Instrument India**

