



Email : sales@laboratoryinstrumentindia.com

Product Name :
CPP White Vinyl Aprons

Product Code :
LBNY-00017-11000470



Description :

Heavy-duty protection in a variety of styles
Highest-quality virgin vinyl in 4, 6, 8, 12 and 18mil thicknesses

Public Utilities:

Cleaning and de-greasing
Industrial maintenance
Janitorial and maintenance work

Benefits:

Superior resistance to chemicals, fats, oils, greases, punctures and abrasions

Optional features:

Stomach patches for extended wear
Raw or sewn edges
Three-piece ties
Rough texture
Food Preparation, Processing and Handling:
Beef/pork/poultry/fish processing
Dairy processing

Chemical:

Chemical processing, oil refining and petrochemicals.

Technical Specification :

CPP White Vinyl Aprons



Laboratory Instrument India