



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Crocodile Clip

**Product Code :**  
LBNY-0003-1010006



**Description :**

Crocodile Clip

**Technical Specification :**

Crocodile Clip mostly use in experimental Instrument cors in labs such as Variable DC Power Supply as well as in many DIY projects for temporary connections of power supply or any other required connection. With its Specific Red and Black insulating cover, one can easily manage to make and Break the connections with no polarity conflicts, also avoids the damage or Short circuits.  
Color: Red and black



**Laboratory Instrument India**