



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microporous Composite Fabric Cleanroom Coveralls, Bulk Packed

Product Code :
LBNY-00017-11000473



Description :

Unique breathable composite material provides particulate and light splash protection for critical environments. Coveralls feature a zipper front, serged seams, and elastic wrists, ankles, and waist.

Made from a microporous composite fabric

Serged seams

Elastic back

Covered elastic wrist and ankles

Bound neck

Lightweight

Anti-static treated

Low linting

Zipper Closure.

Technical Specification :

Microporous Composite Fabric Cleanroom Coveralls, Bulk Packed



Laboratory Instrument India

