



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
120 Series Coveralls

**Product Code :**  
LBNY-00017-11000480



**Description :**

QC is polyethylene-coated fabric provides at least 30 minutes of protection against 42 chemical challenges  
Flexible and durable light liquid splash protection  
Yellow for high visibility  
10 mil thickness  
Laydown collar  
Open wrists  
Open ankles.

**Technical Specification :**

120 Series Coveralls



**Laboratory Instrument India**