



Email : sales@laboratoryinstrumentindia.com

Product Name :
50 Coveralls with Hood

Product Code :
LBNY-00017-11000481



Description :

Made with a microporous film that is laminated to a spund bonded polypropylene substrate
Made from microporous film, 50 garments help protect workers who face on-the-job exposure to a range of non-hazardous aerosols, light liquid splash and dry particles.
Standard hood with elastic around the face
Front zipper closure covered with a storm flap
Serged seams
White
Elastic wrist
Elastic ankles
Thickness: 10 mils.

Technical Specification :

50 Coveralls with Hood



Laboratory Instrument India

